



## Product Overview

The T-Era-20D40P-1550-fiber terahertz photoconductive antenna (THz-PCA) is used to detect wideband terahertz pulses in terahertz time-domain systems. The T-Era-20D40P-1550-fiber THz-PCA is made on high resistive ultra-fast epitaxially grown multi-quantum well InGaAs-InAlAs substrates and is packaged in TeTechS' patent pending terahertz chip fiber coupled enclosure module. The enclosure module houses the THz-PCA with a collimating high-resistive silicon lens attached to the back side of the THz-PCA chip, an FC/APC fiber connector and optical collimating and focusing lenses. The device is packaged in a modular format so that it is easy to change the THz-PCA chip inside the enclosure at a fraction of cost. The device is shipped with the silicon lens aligned and packaged on the back side of the THz-PCA chip. The silicon lens can be re-aligned after changing the THz-PCA chip using our silicon lens setting fixtures.

In the receiving operation mode, the detected THz photocurrent can be measured through the MMCX connector. The standard  $\varnothing 1"$  threaded body makes it convenient to attach the module to other standard optical components or mount it on an optical bench. When excited by optical pulses with 8 mW average optical power, the T-Era-20D40P-1550-fiber THz-PCA generates 200 pA peak terahertz photocurrent with more than 50dB terahertz power spectrum dynamic range.

## Product Specifications

Optical Excitation Wavelength	1540 nm-1560 nm
Average Optical Power	1 mW-15 mW
Bias Voltage for Optical Alignment	5 V
Dark Resistance	1 M $\Omega$
Spectrum Bandwidth	>2.5 THz
Power Spectrum Dynamic Range	>50 dB
Size (O.D., L)	1", 2.5"



## I-V Curves

### I. Dark current versus bias voltage

Figure 1 shows the dark current versus applied bias voltage across the T-Era-20D40P-1550-fiber THz-PCA. A dark resistance around 1 M $\Omega$  is measured.




Figure 1. Dark current versus applied bias voltage across the T-Era-20D40P-1550-fiber THz-PCA.

### Dark current versus bias voltage measurement settings

THz-PCA Under Test	T-Era-20D40P-1550-fiber
Average Optical Power on THz-PCA	0 mW
Bias Voltage on THz-PCA	0.5V-8V

## II. Photocurrent versus bias voltage

Figure 2 shows the photocurrent versus bias voltage across the T-Era-20D40P-1550-fiber THz-PCA. The average optical power on the device is fixed at 20 mW. The photocurrent shows a linear increase with bias voltage.




Figure 2. Photocurrent versus applied bias voltage across the T-Era-20D40P-1550-fiber THz-PCA.

### Photocurrent versus bias voltage measurement settings

THz-PCA Under Test	T-Era-20D40P-1550-fiber
Optical Excitation Wavelength	1550 nm
Optical Pulse Duration	100 fs
Average Optical Power on THz-PCA	8 mW
Bias Voltage on THz-PCA	0.5V-8V

### III. Terahertz Measurement Setup

Figure 3 shows a terahertz response measurement setup for T-Era-20D40P-1550-fiber THz-PCA.




Figure 3. A terahertz response measurement setup for T-Era-20D40P-1550-fiber THz-PCA.

#### Typical THz Time-Domain Measurement Settings

Transmitter Module	T-Era-100A-1550-fiber
Receiver Module	T-Era-20D40P-1550-fiber
Optical Excitation Wavelength	1550 nm
Optical Pulse Duration	100 fs
Average Optical Power on Transmitter	14 mW
Average Optical Power on Receiver	8 mW
Bias Voltage on Transmitter	40 V

#### IV. Terahertz Response

Figure 4 shows a typical THz pulse and its corresponding power spectrum generated by a T-Era-100A-1550-fiber THz-PCAs and detected by a T-Era-20D40P-1550-fiber THz-PCAs in a terahertz time-domain system.




Figure 4. Typical THz pulse and its corresponding power spectrum generated by a T-Era-100A-1550-fiber transmitter module and detected by a T-Era-20D40P-1550-fiber receiver module.

## V. Terahertz peak photocurrent versus bias voltage

Figure 5 shows the terahertz peak photocurrent on the T-Era-20D40P-1550-fiber receiver versus bias voltage on the T-Era-100A-1550-fiber transmitter. The average optical power on the transmitter antenna is 14 mW and on the receiver antenna is 8 mW.




Figure 5. Terahertz peak photocurrent versus bias voltage on the T-Era-20D40P-1550-fiber receiver versus bias voltage on the T-Era-100A-1550-fiber transmitter.